

HELPFUL TIPS FOR HETTICH'S CENTRIFUGES:

The following questions could help you narrow down the correct centrifuge for your application;

~ Tube Type – Application?

Microliter Tubes

Round Bottom Tubes

Conical /Shirted Tubes

Blood Bags

High Capacity Tubes

Plates

~ Tube, Bag or Plate Size?

~ How many per run?

~ At what RCF?

~ Do you need refrigeration?

~ Do you need bio-containment?

~ Do you need a swing out or fixed rotor?

~ Cost / Budget constraints

Note: swing out rotors during centrifugation tubes are swung horizontal

Fixed angle rotors – used primarily in higher speed range

~ RPM – revolutions per minute

